

REFLEX STRAPS

Reflective universal spider straps system for spine boards



Universal retention system for the immobilization of patients compatible with all spinal boards with at least 10 handles

Specific features

- Reflective inserts integrated into the fabric facilitate rescue in conditions of poor visibility
- 4 transversal belts, 3 of which are movable to guarantee adaptation to the spinal board in use
- Width adjustment range thanks to the strap closures

Technical data ¹

Lengths	Min (mm)	Max (mm)
"V" tapes	250	520
Transversal tape 1	750	1350
Transversal tape 2	600	1250
Transversal tape 3	550	1110
Transversal tape 4	480	820
Central belt	1200	
Materials	Nylon	
Weight	630 ± 20 g	

¹ Indicated lengths refer to the installed device. The effective length may vary depending on the device on which the product is applied.

Class I Med.Device according to Dir. 93/42/CEE

REFLEX STRAPS FOR SPINE BOARDS

ST02020A

CND classification

Registration number



V9099

99416

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it