

RES-Q-SPLINT PLUS KIT

Vacuum splints Kit



Vacuum splints for immobilization of limbs. The device allows to limit damage and pain of patients during transfer procedures.

Specific features

PVC sheet easy to clean
 Fast closure thanks to the integrated tear-off straps
 Non-return valve characterized by large air passages allow to reach the desired vacuum level faster
 Carry bag with handles and housing for pump on the rear side
 The supplied pump is characterized by small size and high efficiency

Technical data

Dimensions Size. S	735x640 ± 10 mm
Dimensions Size. M	1100x695 ± 10 mm
Dimensions Size. L	1300x720 ± 10 mm
Bag dimensions	700x500x250 ± 20 mm
Materials	PVC,Al, PS
Weight of the kit	3,8 ± 0,4 kg

Class I Med.Device according to Dir. 93/42/CEE

RES-Q-SPLINT PLUS KIT -
 UNIVERSAL VACUUM SPLINTS

QM22590A

CND Classification



V0804

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it