

RES-Q-SPLINT PLUS

Vacuum splints



Vacuum splints for immobilization of limbs. The device allows to limit damage and pain of patients during transfer procedures.

Specific features

- PVC sheet easy to clean
- Fast closure thanks to the integrated tear-off straps
- Non-return valve characterized by large air passages allow to reach the desired vacuum level faster
- Carry bag with handles and housing for pump on the rear side
- The supplied pump is characterized by small size and high efficiency

Technical data

Dimensions Size. S	735x640 ± 10 mm
Dimensions Size. M	1100x695 ± 10 mm
Dimensions Size. L	1300x720 ± 10 mm
Bag dimensions	700x500x250 ± 20 mm
Materials	PVC,Al, PS
Weight of the kit	3,8 ± 0,4 kg

Class I Med.Device according to Dir. 93/42/CEE

RES-Q-SPLINT PLUS KIT -
UNIVERSAL VACUUM SPLINTS
QM22590A



CND Classification

V0804

Rev.0 (11/05/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it