

RES-Q-SPLINT PLUS

Vacuum splint



Vacuum splint for immobilization of limbs. The device allows to limit damage and pain of patients during transfer procedures.

Specific features

- PVC sheet easy to clean
- Fast closure thanks to the integrated tear-off straps
- Non-return valve characterized by large air passages allow to reach the desired vacuum level faster
- Suitable for leg immobilization

Technical data

Dimensions	1300x720 ± 10 mm
Materials	PVC,Al, PS
Weight	1,05 ± 0,1 kg

Class I Med.Device according to Dir. 93/42/CEE

RES-Q-SPLINT PLUS SIZE L -
UNIVERSAL VACUUM SPL.

QM22591A

CND Classification V0804

Registration number 102931



Rev.2(18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it