

## RES-Q-SPLINT PLUS

### Vacuum splint



Vacuum splint for immobilization of limbs. The device allows to limit damage and pain of patients during transfer procedures.

#### Specific features

- PVC sheet easy to clean
- Fast closure thanks to the integrated tear-off straps
- Non-return valve characterized by large air passages allow to reach the desired vacuum level faster
- Suitable for ankle immobilization

#### Technical data

Dimensions	735x640 ± 10 mm
Materials	PVC,Al, PS
Weight	0,6 ± 0,1 kg

Class I Med.Device according to Dir. 93/42/CEE

RES Q SPLINT PLUS ANKLE  
QM22593A



CND Classification	V0804
Registration number	102927

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)