

RES-Q-SPLINT KIT

Vacuum splints – Kit 3 sizes with pump and bag



Vacuum splints for immobilization of limbs. The device allows to limit damage and pain of patients during transfer procedures.

Specific features

- PVC sheet easy to clean
- Separate channels avoid movement and accumulation of filling material, ensuring its better distribution making more effective the immobilization of the limb.
- Fast closure thanks to the integrated tear-off straps
- Non-return valve characterized by large air passages allow to reach the desired vacuum level faster
- Carry bag with handles and housing for pump on the rear side
- The supplied pump is characterized by small size and high efficiency

Technical data

Dimensions Size. S	560x460 ± 10 mm
Dimensions Size. M	490x700 ± 10 mm
Dimensions Size. L	750x900 ± 10 mm
Bag dimensions	550x410x160 ± 20 mm
Materials	PVC,Al, PS
Weight of the kit	3,16 ± 0,2 kg

Class I Med.Device according to Dir. 93/42/CEE

RES-Q-SPLINT - KIT 3 SIZES
W/PUMP AND BAG

QM22500A

CND Classification V0804

Registration number 195244



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it