

# RES-Q-SPLINT

## Vacuum splints



Vacuum splints for immobilization of limbs. The device allows to limit damage and pain of patients during transfer procedures.

### Specific features

- PVC sheet easy to clean
- Separate channels avoid movement and accumulation of filling material, ensuring its better distribution making more effective the immobilization of the limb.
- Fast closure thanks to the integrated tear-off straps
- Non-return valve characterized by large air passages allow to reach the desired vacuum level faster
- Carry bag with handles and housing for pump on the rear side
- The supplied pump is characterized by small size and high efficiency

### Technical data

Dimensions Size. S	560x460 ± 10 mm
Dimensions Size. M	490x700 ± 10 mm
Dimensions Size. L	750x900 ± 10 mm
Bag dimensions	550x410x160 ± 20 mm
Materials	PVC,Al, PS
Weight of the kit	3,16 ± 0,2 kg

Class I Med.Device according to Dir. 93/42/CEE

RES-Q-SPLINT - KIT 3 SIZES  
W/PUMP AND BAG  
QM22500A



CND Classification	V0804
Registration number	195244

Rev.0 (11/05/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)