

RES-Q-SPLINT

Vacuum splint



Vacuum splints for immobilization of limbs. The device allows to limit damage and pain of patients during transfer procedures.

Specific features

- PVC sheet easy to clean
- 3 separate channels, limit displacement and accumulations of the filling material, ensuring a better distribution and making the immobilization more effective
- Fast closure thanks to the integrated tear-off straps
- Non-return valve characterized by large air passages allow to reach the desired vacuum level faster

Technical data

Longer side	490 ± 10 mm
Shorter side	380 ± 10 mm
Lenght	700 ± 10 mm
Material	PVC,Al, PS
Weight	630 ± 50g

Class I Med.Device according to Dir. 93/42/CEE

RES-Q-SPLINT SIZE M - VAC.

SPL.WITH SEPAR. CHAMB.

QM22520A

CND Classification V0804

Registration number 102924



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it