

## RES-Q-SPLINT

### Vacuum splint



Vacuum splints for immobilization of limbs. The device allows to limit damage and pain of patients during transfer procedures.

#### Specific features

- PVC sheet easy to clean
- Fast closure thanks to the integrated tear-off straps
- Non-return valve characterized by large air passages allow to reach the desired vacuum level faster

#### Technical data

Dimensions Size. S	560x460 ± 10 mm
Materials	PVC,Al, PS
Weight	295 ± 50g

Class I Med.Device according to Dir. 93/42/CEE

RES-Q-SPLINT SIZE S - VAC.SPL.  
WITH SEPAR. CHAMB.

QM22530A



CND Classification V0804

Registration number 102933

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)