

EMSoft

Vacuum splints



Immobilization splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations. Ideal for wrist immobilization.

Specific Characteristics

- Low storage space required
- One-way valve
- Vinyl layer to make cleaning procedures easier
- High moldability
- Fast application thanks to the strap closures
- Vacuum pump and carry bag included

Technical data

	Dimensions (mm)	Weight (g)
Size S	750x600	585
Size M	900x700	830
Size L	1100x800	1130
Pump	Ø36 x 310	500

Note: Dimensions tolerance ± 10mm, Weight tolerance ± 40g

Class I Med.Device according to Dir. 93/42/CEE

EMSOFT KIT THREE SIZES -
VACUUM SPLINTS
QM22707A



CND Classification	V0804
Registration number	195248

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it