

# EMSoft

## Vacuum splints



Immobilization splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations. Ideal for wrist immobilization.

### Specific Characteristics

- Low storage space required
- One-way valve
- Vinyl layer to make cleaning procedures easier
- High moldability
- Fast application thanks to the strap closures
- Vacuum pump and carry bag included

### Technical data

	Dimensions (mm)	Weight (g)
Size S	750x600	585
Size M	900x700	830
Size L	1100x800	1130
Pump	Ø36 x 310	500

Note: Dimensions tolerance ± 10mm, Weight tolerance ± 40g

Class I Med.Device according to Dir. 93/42/CEE

EMSOFT KIT THREE SIZES -  
VACUUM SPLINTS  
QM22707A



CND Classification	V0804
Registration number	195248

Rev.1 (30/06/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)