

ECO AIR SPLINT

Inflatable Splints



Inflatable splints designed for immobilization. They allow to limit patient pain and damages that may occur during rescue operations.

Specific features

- Low storage space required
- Quick application thanks to the zip closure
- Varying the inflation pressure is possible to increase or decrease the immobilization level of the limb
- Color coded zip puller for quick size identification

Technical data

Included sizes	Arm 37 cm
	Leg 40 cm
	Arm 64 cm
	Leg 67 cm
	Arm 75 cm
	Leg 85 cm
Material	PVC
Weight	1500 ± 50 g

Class I Med.Device according to Dir. 93/42/CEE

ECO AIR SPLINT - KIT 6 SIZES WITH BAG

QS24200A

CND Classification



V0804

Registration number

195271

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it