

BLU SPLINT

Rigid splint with flexible internal structure



Immobilization splint for traumatized limbs. Used to limit the damage and pain to the patient related to rescue operations. Ideal for leg immobilization

Specific features

- Flexible and moldable core
- Strap closure
- Washable up to 40 ° C by removing the inner core
- Soft materials to ensure high comfort

Technical data

Splint	Size ⁽¹⁾ (mm)	Weight (g)
Leg	580x450	360
Materials	Neoprene, Al, Nylon	

⁽¹⁾ Measures taken with expanded splint and expressed as Length x maximum width

Class I Med.Device according to Dir. 93/42/CEE

BLUE SPLINT LEG

JM80003A



CND Classification

V0804

Rev.1 (17/01/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it