

BLUE SPLINT

Rigid splint with flexible internal structure



Blue Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- Flexible and moldable core
- Strap closure
- Washable up to 40 °C by removing the inner core
- Soft materials to ensure high comfort

Technical data

Splint	Size ⁽¹⁾ (mm)	Weight (g)
Wrist	310 x 240 ± 10	145
Materials	Neoprene, Al, Nylon	

⁽¹⁾ Measures taken with expanded splint and expressed as Length x maximum width

Class I MD compliant with UE Reg. 2017/745

BLUE SPLINT – RIGID SPLINT FOR WRIST

JM80006A

CND Classification



V0804

Registration number

103035

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it