

BLU SPLINT PRO

Rigid splint with flexible internal structure



Immobilization splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations. Ideal for forearm immobilization.

Specific Characteristics

- Flexible and moldable core
- Quick release fastener
- Suitable to be used as an armrest during intravenous infusions
- Washable up to 40 ° C by removing the inner core
- Soft materials to ensure high comfort

Technical data

| | |
|------------|---------------------|
| Dimensions | 390x300 mm |
| Materials | Neoprene, Al, Nylon |
| Weight | 270 g |

⁽¹⁾ Measures taken with expanded splint and expressed as Length x maximum width

Class I Med.Device according to Dir. 93/42/CEE

BLUE SPLINT PRO FOREARM -
RIGID SPLINTS WITH FLEXIBLE
INTERNAL STRUCTURE
JM80032A



| | |
|---------------------|--------|
| CND Classification | V0804 |
| Registration number | 103045 |

Rev.0 (08/06/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it