

BLU SPLINT PRO

Rigid splint with flexible internal structure



Immobilization splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations. Ideal for elbow and ankle immobilization.

Specific characteristics

- Flexible and moldable core
- Quick release fastener
- Suitable to be used as an armrest during intravenous infusions
- Washable up to 40 ° C by removing the inner core
- Soft materials to ensure high comfort

Technical data

Dimensions	540x280 mm
Materials	Neoprene, Al, Nylon
Weight	390 g

⁽¹⁾ Measures taken with expanded splint and expressed as Length x maximum width

Class I Med.Device according to Dir. 93/42/CEE

BLUE SPLINT PRO KIT - RIGID
SPLINTS WITH FLEXIBLE INTERNAL
STRUCTURE
JM80034A



CND Classification	V0804
Registration number	103053

Rev.0 (08/06/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it