

# FIXO SPLINT

Rigid splint



Rigid splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations.

### Specific Characteristics

- Low storage space required
- Quick application thanks to the strap belt system
- Heat sealed coating to ensure easy sanitation
- Soft padding to increase patient comfort

### Technical data

Dimensions <sup>(1)</sup>	150x230x730 mm
Materials	PVC, wood
Weight	870 g

<sup>(1)</sup>Width x Depth x Height (as shown in the picture)

Class I Med.Device according to Dir. 93/42/CEE

FIXO SPLINT SIZE L - RIGID SPLINTS

QC70050A



CND Classification V0804

Registration number 103657

Rev.0 (08/06/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)