

FIXO SPLINT

Rigid splint



Rigid splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations.

Specific Characteristics

- Low storage space required
- Quick application thanks to the strap belt system
- Heat sealed coating to ensure easy sanitation
- Soft padding to increase patient comfort

Technical data

Dimensions ⁽¹⁾	150x230x730 mm
Materials	PVC, wood
Weight	870 g

⁽¹⁾ Width x Depth x Height (as shown in the picture)

Class I Med.Device according to Dir. 93/42/CEE

FIXO SPLINT SIZE L - RIGID SPLINTS

QC70050A



CND Classification	V0804
Registration number	103657

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it