

FIXO SPLINT KIT

Rigid splints



Rigid splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations .

Specific Characteristics

- Low storage space required
- Quick application thanks to the strap belt system
- Heat sealed coating to ensure easy sanitation
- Soft padding to increase patient comfort

Standard equipment

- | | |
|----------|----------------------|
| QC70100A | Orange PVC carry bag |
|----------|----------------------|

Technical data

Splint	Dimensions ⁽¹⁾ (mm)	Weight (g)
Size S	120x230x500	640
Size L	150x230x730	870
Materials	PVC, wood	
Weight of the kit	2 kg	
Bag dimensions	800x430 mm	
Bag weight	450 g	

Class I Med.Device according to Dir. 93/42/CEE

FIXO SPLINT KIT TWO SIZES - RIGID
SPLINTS

QC70090A

CND Classification V0804

Registration number 195300



⁽¹⁾ Width x Depth x Height (as shown in the picture)

Rev.1 (19/11/2014)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it