

# NEXUS

## Anti-shrinkage vacuum mattress



Vacuum mattresses are devices for the immobilization of trauma patients. The vacuum system allows different levels of immobilization. The more rigid helps the right positioning of the patient; softer ones allow to mitigate the low-frequency vibrations typically found on emergency vehicles during transport.

### Specific characteristics

Robust, non-absorbent and self-extinguishing fabric  
 Reinforced handles to ensure a safe grab  
 Embedded inserts in the lower side of the mattress, offer more support to the patient and limit the retraction of the mattress during decompression procedures

### Standard equipment

STX 597 - ONE PIECE PLASTIC BELT

### Accessories

QM22109A ABS double effect pump  
 QM22120A Aluminum Pump  
 QM22170A Carry bag for vacuum mattresses  
 QM22199A Repair kit

### Technical data

Length	2025 ± 10mm
Width	855 ± 10mm
Thickness (valve zone)	100 ± 10mm
N° of handles	8
Materials	PVC, Wood, Nylon Al
Weight <sup>With belts</sup>	5,6 ± 0,5 kg
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

NEXUS - ANTI-SHRINKAGE  
 VACUUM MATTRESS  
**QM22900A**



**CND Classification** V0804  
**Registration number** 102669

Rev.2 (18/07/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)