

# NEXUS

## Anti-shrinkage vacuum mattress



Vacuum mattresses are devices for the immobilization of trauma patients. The vacuum system allows different levels of immobilization. The more rigid helps the right positioning of the patient; softer ones allow to mitigate the low-frequency vibrations typically found on emergency vehicles during transport.

### Specific characteristics

Robust, non-absorbent and self-extinguishing fabric  
Reinforced handles to ensure a safe grab  
Embedded inserts in the lower side of the mattress, offer more support to the patient and limit the retraction of the mattress during decompression procedures

### Standard equipment

ST00597A	STX 597 - ONE PIECE PLASTIC BELT
QM22120A	VACUUM PUMP MOD. 120
QM22902A	TRANSPORT BAG

### Accessories

QM22109A	ABS double effect pump
QM22120A	Aluminum Pump
QM22170A	Carry bag for vacuum mattresses
QM22199A	Repair kit

### Technical data

Length	2025 ± 10mm
Width	855 ± 10mm
Thickness (valve zone)	100 ± 10mm
N° of handles	8
Materials	PVC, Wood, Nylon Al
Weight <sup>With belts</sup>	5,6 ± 0,5 kg
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

NEXUS - ANTI-SHRINKAGE  
VACUUM MATTRESS



**QM22901A**

<b>CND Classification</b>	V0804
<b>Registration number</b>	102670

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)