

VI MATT 2

Veterinary vacuum mattress



Veterinary vacuum mattress intended to immobilize small sized pets. The device can also be used to provide greater comfort during dressing procedures.

Specific characteristics

Robust, non absorbent and easy to clean fabric
 Extremely low weight
 Valve with safety system to limit loss of vacuum during application of the device

Technical data

Length	600 ± 30 mm
Width	300 ± 30 mm
Materials	PVC, PS
Weight	240 ± 30 g

Accessories

QM22109A	QMX 109 - hand vacuum pump abs
QM22120A	QMX 120 - hand vacuum pump aluminium
QM22125A	QMX 125 - hand vacuum pump aluminium
QM22130A	QMX 130 - hand vacuum pump al. w/logo

SMALL VACUUM MATTRESS
 CM.30X60 IN BLUE PVC
QM22411A

Rev.0 (14/07/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it