

## VI MATT 1

### Veterinary vacuum mattress



Veterinary vacuum mattress intended to immobilize medium sized pets. The device can also be used to provide greater comfort during dressing procedures.

#### Specific characteristics

Robust, non absorbent and easy to clean fabric  
Extremely low weight  
Valve with safety system to limit loss of vacuum during application of the device

#### Technical data

Length	1200 ± 30 mm
Width	600 ± 30 mm
Materials	PVC, PS
Weight	0,85 ± 0,3 kg

#### Accessories

QM22109A	QMX 109 - hand vacuum pump abs
QM22120A	QMX 120 - hand vacuum pump aluminium
QM22125A	QMX 125 - hand vacuum pump aluminium
QM22130A	QMX 130 - hand vacuum pump al. w/logo

VACUUM MATTRESS CM. 60X120  
IN BLUE PVC  
**QM22412A**

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)