

SPENCER CONTOUR

Anatomic universal head immobilizer



Device for head immobilization intended to be used with Spencer spinal boards or pickup stretchers, except pediatric models.

Specific features

- Openings for ear inspection
- Soft pillows to increase comfort
- Wide strap surfaces make the headblock positioning easier
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Smooth surfaces for easy sanitation
- Central removable padded area for adult or pediatric use

Technical data

Pillow dimensions ¹	235x130x160 mm
Base dimensions ²	390x260x15 mm
Head immobilizer weight	1,5 ± 0,1 kg
Materials	PU, Nylon, PVC
¹ One pillow dimensions	
² Buttonhole excluded	

Reference standard

NA

Class I Med.Device according to Dir. 93/42/CEE

SPENCER CONTOUR BLUE - ANAT.

UNIV. HEAD IMMOB.

SH00200A

CND Classification



Y060312

Rev.2 (25/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it