

# SPENCER CONTOUR

Anatomic universal head immobilizer



Device for head immobilization intended to be used with Spencer spinal boards or pickup stretchers, except pediatric models.

### Specific features

- Openings for ear inspection
- Soft pillows to increase comfort
- Wide strap surfaces make the headblock positioning easier
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Smooth surfaces for easy sanitation
- Central removable padded area for adult or pediatric use

### Technical data

Pillow dimensions <sup>1</sup>	235x130x160 mm
Base dimensions <sup>2</sup>	390x260x15 mm
Head immobilizer weight	1,5 ± 0,1 kg
Materials	PU, Nylon, PVC
<sup>1</sup> One pillow dimensions	
<sup>2</sup> Buttonhole excluded	

### Reference standard

NA

Class I Med.Device according to Dir. 93/42/CEE

SPENCER CONTOUR  
YELLOW/BLACK - AN.UN.HEAD  
IMM.



SH00201A

CND Classification Y060312

Registration number 102237

Rev.2 (30/06/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)