

SUPER BLACK

Compact universal head immobilizer



Device for head immobilization intended to be used with Spencer spinal boards or pickup stretchers. It can be used along with other immobilization devices like SED and all Spencer cervical collars.

Specific Characteristics

- The padding of the base increase patient comfort
- Wide strap surfaces make the headblock positioning easier
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Pillows made of non-absorbent material for an easy sanitation
- Large holes for ear inspection

Technical data

Base dimensions ^(buttonhole excluded)	365x260x20 mm
Pillow dimensions ^(1 piece)	255x85x160 mm
Hole for ear inspection	ø70 mm
Base weight	360 ± 10g
Pillow weight	500 ± 10g
Overall weight	1,46 ± 0,07 kg
Materials	PE, PVC, Nylon

Class I Med.Device according to Dir. 93/42/CEE

SUPER BLUE BLACK - COMPACT

UNIV. HEAD IMMOB.

SH00104C

CND Classification Y060312

Registration number 161552



Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it