

## SUPER BLUE HP

universal head immobilizer



Device for head immobilization intended to be used with Spencer spinal boards or pickup stretchers. It can be used along with other immobilization devices.

### Specific Features

- Wide strap surfaces make the headblock positioning easier
- On the lower side of the base, a large strap element covers the entire surface allowing a more stable fixing on the immobilization device
- Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
- Pillows made of non-absorbent material for an easy sanitation
- Large holes for ear inspection

### Technical data

Base dimensions	445x280x15 mm
Pillow dimensions <sup>(1 piece)</sup>	255x80x155 mm
Hole for ear inspection	ø70 mm
Base weight	700 ± 10g
One pillow weight	500 ± 30g
Overall weight	1,83 ± 0,3 kg
Materials	PE, PVC, Nylon

Class I Med.Device according to Dir. 93/42/CEE

MOD.F011 SUPER BLUE HP  
UNIVERSAL HEAD IMMOBILIZER  
SH00166C



CND Classification Y060312

Rev.0 (26/07/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)