

FXA PRO

Anatomic universal ankle immobilizer



Device for ankle immobilization used with all Spencer spinal boards.

Specific Characteristics

The wrap-around shape of the support area makes it imperceptible to the patient

One-piece adjustable strap band for device fixing and for ankle immobilization

The band can be removed allowing easy cleaning procedures

Extremely robust and lightweight

Technical data

Dimensions	290x160x80 ± 5mm
Materials	PE, Nylon
Weight	350 ± 10g

Class I Med.Device according to Dir. 93/42/CEE

FXA PRO - ANATOMIC UNIVERSAL
ANKLE IMMOBILIZER

SH00151A

CND Classification V0899

Registration number 195166



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it