

## QMX 109

Pump for immobilization devices



Double effect pump designed for use with Spencer vacuum mattresses or vacuum/inflatable splints.

### Specific features

- Large air chamber allow faster depression/inflation
- The shape of the base allow to hold firmly the pump with the feet during use
- Handles with non-slip finish
- The suction/inflation tube is connected through a quick lock/release system
- Connections with clear identification on the handle group (inflate/deflate)

### Technical data

Dimensions	215x100x445 ± 5 mm
Tube length	100 ± 5 cm
Materials	ABS, PE
Weight	0,6 ± 0,1 kg

HAND VACUUM PUMP ABS  
 SPENCER MOD. 109  
**QM22109A**

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)