

QMX 120

Pump for vacuum systems



Pump intended for use with with Spencer vacuum mattresses or splints.

Specific features

The aluminum chamber makes it a rugged and lightweight device

The shape of the base allows to hold the pump firmly during use

Technical data

Dimensions	500x140x80 mm
Tube length	100 ± 5 cm
Materials	Al, PVC, Steel
Weight	1,14 ± 0,1 kg

ALUMINUM VACUUM PUMP
QMX120
QM22120A

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it