

## QMX 130

### Aluminum vacuum pump



Compact vacuum pump, ideal for use with vacuum splints.

#### Specific features

The aluminum chamber makes the device strong and light  
 Characterized by small dimensions, it can be easily stored  
 in any vacuum splints kit  
 Rounded knob for firm grip

#### Technical data

Dimensions	Ø55x320 mm
Tube internal diameter	8 mm
Tube length	255 ± 10mm
Materials	Al, PVC
Weight	400 ± 10 g

QMX 130 - ALUMINIUM PUMP FOR  
 VACUUM SPLINTS  
**QM22130A**

Rev.1 (29/08/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)