

## QMX 199

### Repair kit



Repair kit for vacuum mattresses or vacuum splints.

#### Specific Characteristics

4 PVC patches  
1 tube of glue

#### Technical data

One patch dimensions	ø5 mm
Glue	5 g
Overall weight	14 g

REPAIR KIT  
QMX22199A

Rev.0 (06/10/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)