

QMX 199

Repair kit



Repair kit for vacuum mattresses or vacuum splints.

Specific Characteristics

4 PVC patches
1 tube of glue

Technical data

One patch dimensions	ø5 mm
Glue	5 g
Overall weight	14 g

REPAIR KIT
QMX22199A

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it