

CROSS

10G Self loading stretcher



Spencer Cross is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. The loading carriage can be folded to make movement in confined spaces easier. 10G certified when used with the fastener Spencer S-Max.

Caratteristiche specifiche

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- Trendelenburg antishock position
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bars *Safe bar*
- Adjustable heights system to facilitate patient transfer procedures
- Foldable loading carriage to make handling in confined spaces easier

Optional accessories

| | |
|----------|--|
| ST42707C | S-Max fastener <small>(Mandatory for 10g versions)</small> |
| ST70000A | PVC anatomic mattress QMX/777 |
| EN90003C | END-T Table for stretchers 10g cert. |
| ST00491A | STX 90 Telescopic headrest for stretchers |
| ST00497B | DNA STRAP Belt w/integrated re-winding system |
| ST00498B | DNA STRAP Thor.Belt w/integr. re-winding syst. |
| IF01047C | Track 4-30 telescopic IV holder |
| CB09025C | Tanker oxygen tank holder for stretcher side bar |
| CB09028C | Tanker support for Safe Bar |

Standard equipment

| | |
|----------|---|
| ST70002A | STX 702 two pieces black reflex metal belts |
|----------|---|

Technical data

| | |
|-------------------------------|------------------|
| Length | 1980 ± 10 mm |
| Shortened length | 1660 ± 10 mm |
| Width | 580 ± 10 mm |
| Weight | 45 ± 1 kg |
| Loading capacity | 170 kg |
| Wheels | ø200 mm ± 5% |
| Height loading carriage wheel | From 52 to 70 cm |
| Frame material | Steel |
| Side bars length | 680 mm |
| Side bars Height | 200 mm |

Class I Med.Device according to Dir. 93/42/CEE

CROSS T SELF LOADING STRETCHER
10G YELLOW BOARD BLACK
FRAME



CR60045B

| | |
|----------------------------|--------|
| CND Classification | V0899 |
| Registration number | 151557 |

Rev.0 (22/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it