

CROSSOVER

10g Self Loading stretcher



Spencer Crossover is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. The loading carriage can be folded to make movement in confined spaces easier. Crossover is characterized by higher loading capacity than a standard Cross.

Specific Characteristics

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- Antishock positions Trendelenburg
- Extreme maneuverability thanks to the Twist system which, when operated, makes pivoting also the front wheels.
- Brakes on rear wheels
- Foldable side bars Safe bar
- Adjustable heights system
- Foldable loading carriage

Accessories

| Accessories | |
|-------------|---|
| ST42702C | R-Max fixing system <small>(Mandatory for 10g versions)</small> |
| EN90003C | END-T Table for stretchers 10g cert. |
| ST00491A | STX 90 Telescopic headrest for stretchers |
| ST00497B | DNA STRAP Belt w/integrated re-winding system |
| ST00498B | DNA STRAP Thorax. Belt w/integrated re-winding system |
| IF01047C | Track 4-30 telescopic IV holder |
| CB09025C | Tanker oxygen tank holder for stretcher side bar |
| CB09028C | Tanker support for Safe Bar |

Standard Equipment

- PVC anatomic mattress QMX/777
- STX 702 two pieces black reflex metal belts

Technical data

| | |
|-------------------------------|------------------|
| Length | 1980 mm |
| Shortened length | 1660 mm |
| Width | 580 mm |
| Weight | 49 kg |
| Loading capacity | 250 kg |
| Wheels | ø200 mm |
| Height loading carriage wheel | From 52 to 70 cm |
| Frame material | Stainless Steel |
| Side bars length | 680 mm |
| Side bars Height | 200 mm |

Class I Med.Device according to Dir. 93/42/CEE

CROSS OVER T CHROME FRAME
YELLOW BOARD W/QMX 777 10G

CR70004C

| | |
|---------------------|--------|
| CND Classification | V0899 |
| Registration number | 380192 |



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it