

# CROSSOVER

## 10g Self Loading stretcher



Spencer Crossover is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. The loading carriage can be folded to make movement in confined spaces easier. Crossover is characterized by higher loading capacity than a standard Cross. 10G certified when used with the fastener Spencer S-Max

### Specific Characteristics

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- Antishock positions Trendelenburg
- Extreme maneuverability thanks to the Twist system which, when operated, makes pivoting all 4 wheels
- Brakes on rear wheels
- Foldable side bars Safe bar
- Adjustable heights system
- Foldable loading carriage

### Standard Equipment

- PVC anatomic mattress QMX/777
- STX 702 two pieces black reflex metal belts

### Technical data

Length	1980 ± 10 mm
Shortened length	1660 ± 10mm
Width	570 ± 10mm
Weight	45 kg
Loading capacity	250 kg
Wheels	ø200 mm
Height loading carriage wheel	From 47 to 55 cm
Frame material	Stainless steel
Side bars length	680 mm
Side bars Height	200 mm

### Accessories

ST42707C	S-MAX 10g fastener for stretchers <small>(Mandatory for 10g versions)</small>
ST00499B	STX 499 Thorax. Belt
IF01047C	Track 4-30 telescopic IV holder
ST00479C	FIXO KID intended for secure fixation of various age children
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thor.Belt w/integr. re-winding syst.
ST00491A	STX 90 Telescopic headrest for stretchers

### Class I Med.Device according to Dir. 93/42/CEE

LOW CROSSOVER T F/CHROME  
P/BLACK W/QMX777 10G CER.



### CND Classification

V0899

Rev.0 (26/04/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)