


CROSSOVER SHARP

10G Self-loading stretcher with integrated quick release system



Spencer Crossover Sharp is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. The loading carriage can be folded to make movement in confined spaces easier. Crossover Sharp integrates a manual system strategically placed on the commands handlebar that allows to unlock the stretcher from the fastener S-MAX.

Specific features

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- Trendelenburg antishock position
- 2 rear large diameter pivoting wheels to increase sliding action
- Extreme maneuverability thanks to the Twist system which, when operated, makes pivoting also the front wheels.
- Brakes on rear wheels
- Foldable side bars *Safe bar*
- Adjustable heights system makes the transfer of the patient easier
- Foldable loading carriage allows to reduce maneuvering space
- Thanks to the integration of the Sharp system, is no longer necessary to operate on the fastener to unlock the stretcher
- Note: The Sharp system is compatible only with S-MAX fastener

Optionals

ST42707B	10g Fastener S-Max
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integr. re-winding system
IF01047C	Track 4-30 telescopic IV holder
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Standard equipment

- PVC anatomic mattress QMX/777
- STX 702 two pieces black reflex metal belts
- Sharp System

Technical data

Length	1980 mm
Shortened length	1660 mm
Width	580 mm
Weight without standard equipment	45 kg
Loading capacity	250 kg
Wheels	ø200 mm
Height loading carriage wheel	From 52 to 70 cm
Frame material	Stainless Steel
Side bars length	680 mm
Side bars Height	200 mm

Class I Med.Device according to Dir. 93/42/CEE

CROSSOVER SHARP T T/CHROME
B/BLACK W/QMX777 CERT.10G

CR70014



CND Classification

V0899

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza(PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it