

## CROSSOVER

### 10g Self Loading stretcher



Spencer Crossover is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. The loading carriage can be folded to make movement in confined spaces easier. Crossover is characterized by higher loading capacity than a standard Cross.

#### Specific Characteristics

Polyethylene board easy to clean
Color coded commands
Adjustable backrest
Antishock positions Trendelenburg
Extreme maneuverability thanks to the Twist system which, when operated, makes pivoting also the front wheels.
Brakes on rear wheels
Foldable side bars Safe bar
Adjustable heights system
Foldable loading carriage

#### Accessories

ST42702C	R-Max fixing system <small>(Obligatory for 10g versions)</small>
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integrated re-winding system
IF01047C	Track 4-30 telescopic IV holder
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

#### Standard Equipment

ST70000A	PVC anatomic mattress QMX/777
ST70002A	STX 702 two pieces black reflex metal belts

#### Technical data

Length	1980 mm
Shortened length	1660 mm
Width	580 mm
Weight <small>(Without standard equipment)</small>	45 ± 1 kg
Loading capacity	250 kg
Wheels	ø200 mm
Height loading carriage wheel	From 52 to 70 cm
Frame material	Stainless steel
Side bars length	680 mm
Side bars Height	200 mm

#### Class I Med.Device according to Dir. 93/42/CEE

CROSS OVER T CHROME FRAME  
BLACK BOARD WITH QMX777 10G

**CR70002C**

<b>CND Classification</b>	V0899
<b>Registration number</b>	343886



Rev.5 (14/07/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)