

## CARRERA PRO

Self-Loading stretcher

With housing for spine board



Carrera Pro is a self-loading stretcher for patient transport usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands and they open automatically during the unloading procedure. The spinal board can be placed in a dedicated housing located under the patient board.

### Specific Characteristics

Polyethylene board easy to clean
Color coded commands
Adjustable backrest
2 fixed wheels and 2 pivoting with large diameter to make movement easier
Brakes on rear wheels
Foldable side bars <i>Safe bar</i>
Spinal board easy to remove, making immobilization procedures faster

### Accessories

ST42107A	FP Rear fastener for Carrera Pro
ST42710A	FA HUX Anterior fastener for stretchers
ST42200A	FA Anterior fastener for self loading stretchers
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integrated re-winding system
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar
ST70018A	PVC anatomic mattress QMX/776
ST02061B	B-BAK PIN YELLOW SPINE BOARD W/PINS
<b>Note</b>	Is possible to use only the spinal board Spencer B-Bak PIN not included

### Standard Equipment

ST70002A	STX 702 two pieces black reflex metal belts
----------	---

### Technical data

Length	1980 ± 10 mm
Width	585 ± 10 mm
Weight	36 ± 1 kg
Loading capacity	170 kg
Wheels	∅200 mm
Height loading carriage wheel	61 cm
Frame material	Steel
Side bars length	680 mm
Side bars Height	200 mm

### Class I Med.Device according to Dir. 93/42/CEE

CARRERA S PRO HIGH LOADING  
CHROME FRAME BLACK BOARD  
CA01051E



<b>CND Classification</b>	V0899
<b>Registration number</b>	148434

Rev.0 (04/04/2019)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)