

CARRERA S TEC HIGH

10g Self Loading stretcher



Spencer Carrera S Tec is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. 10G certified when used with the fixing system Spencer R-Max or S-MAX.

Specific Characteristics

Polyethylene board easy to clean
Color coded commands
Adjustable backrest
2 fixed wheels and 2 pivoting with large diameter to make movement easier
Brakes on rear wheels
Foldable side bars <i>Safe bar</i>

Accessories

ST42702C	R-Max fixing system <small>(Mandatory for 10g versions)</small>
ST42707C	S-Max fixing system <small>(Mandatory for 10g versions)</small>
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integrated re-winding system
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Standard Equipment

ST70002A	STX 702 two pieces black reflex metal belts
ST70018A	PVC anatomic mattress QMX/776

Technical data

Length	1980 ± 10 mm
Width	580 ± 10 mm
Weight	38 ± 1 kg
Loading capacity	170 kg
Wheels	ø200 mm
Height loading carriage wheel	60 cm
Frame material	Steel
Side bars length	680 mm
Side bars Height	200 mm

Class I Med.Device according to Dir. 93/42/CEE

CARRERA S TEC HIGH LOADING
YELLOW/F BLACK/B WITH
MATTRESS CERTIFIED 10 G



CA71054C

CND Classification	V0899
Registration number	343885

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it