

CARRERA T ACTIV

Self Loading stretcher



Spencer Carrera Activ is a self-loading stretcher for transport of patients. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The front part of the frame can be folded, reducing considerably the length of the stretcher when it is needed. The legs fold away independently, activating the dedicated commands during loading procedures, and open automatically when unloading.

Specific Characteristics

| |
|---|
| Polyethylene board easy to clean |
| Color coded commands |
| Adjustable backrest |
| Board with Trendelenburg antishock position |
| 2 fixed wheels and 2 pivoting with large diameter to make movement easier |
| Brakes on rear wheels |
| Foldable side bars <i>Safe bar</i> |
| Yellow painted frame |
| Can be shortened in less than 2 seconds |
| Greater maneuverability in tight spaces than a standard stretcher |

Standard Equipment

| | |
|----------|---|
| ST70002A | STX 702 two pieces black reflex metal belts |
|----------|---|

Technical data

| | |
|-------------------------------|-----------------|
| Length | 1975 ± 10 mm |
| Width | 580 ± 10 mm |
| Length shortened | 1680 ± 20 mm |
| Weight | 35 kg |
| Loading capacity | 170 kg |
| Wheels | ø200 mm |
| Height loading carriage wheel | 60 cm |
| Frame material | Stainless steel |
| Side bars length | 680 mm |
| Side bars Height | 200 mm |

Accessories

| | |
|----------|---|
| ST42100A | Posterior fastener for self loading stretchers |
| ST42200A | Anterior fastener for self loading stretchers |
| EN90003C | END-T Table for stretchers 10g cert. |
| ST00491A | STX 90 Telescopic headrest for stretchers |
| ST00497B | DNA STRAP Belt w/integrated re-winding system |
| ST00498B | DNA STRAP Thorax. Belt w/integrated re-winding system |
| IF01047C | Track 4-30 telescopic IV holder |
| IF01049B | Track 5 telescopic IV holder for Carrera |
| ST70000A | QMX777 Black PVC anatomic mattress |
| CB09025C | Tanker oxygen tank holder for stretcher side bar |
| CB09028C | Tanker support for Safe Bar |

Class I Med.Device according to Dir. 93/42/CEE

CARRERA T ACTIV HIGH
AUTOLOADING STRETCHER
YELLOW FRAME BLACK BOARD
CA21154C



| | |
|----------------------------|--------|
| CND Classification | V0899 |
| Registration number | 148417 |

Rev.0 (23/03/15)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it