

# CARRERA T TEC

## Self-loading stretcher



Spencer Carrera T Tec is a compact and lightweight self-loading stretcher for transport of patients. It can be used for rescue operations on road and sanitary vehicles. The legs fold independently operating the specific commands and open automatically when unloading.

### Specific Characteristics

Polyethylene board easy to clean
Color coded commands
Adjustable backrest
Patient board with Trendelenburg antishock position
2 fixed wheels and 2 pivoting with large diameter to make movement easier
Brakes on rear wheels
Foldable side bars <i>Safe bar</i>
Yellow painted frame

### Accessories

ST42712A	HUX INTEGRAL FASTENER
ST42100A	FP Posterior fastener for self loading stretchers
ST42200A	FA Anterior fastener for self loading stretchers
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integ. re-winding syst.
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
ST70000A	QMX777 Black PVC anatomic mattress
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

### Standard Equipment

STX 702 two pieces black reflex metal belts
QMX776 PVC anatomic mattress

### Technical data

Length	1980 ± 10 mm
Width	580 ± 10 mm
Weight	37 ± 1kg
Loading capacity	170 kg
Wheels	ø200 mm
Height of the loading carriage wheel	50 cm
Frame material	Steel
Side bars length	680 mm
Side bars Height	200 mm

### Class I Med.Device according to Dir. 93/42/CEE

CARRERA T TEC LOW YELLOW  
FRAME BLACK BOARD MATTRESS  
QMX776



**CA50154E**

<b>CND Classification</b>	V0899
<b>Registration number</b>	180074

Rev.0 (06/02/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)