

ELLA SELF

Self loading stretcher



Spencer Ella Self is a compact and lightweight self-loading stretcher for transport of patients. It can be used for rescue operations on road and sanitary vehicles. The legs fold independently operating the specific commands and open automatically when unloading.

Specific Characteristics

Polyethylene board easy to clean
 Color coded commands
 Adjustable backrest in 5 positions in addition to the horizontal
 Board with Trendelenburg antishock position
 2 fixed wheels and 2 pivoting with large diameter to make movement easier
 Brakes on rear wheels
 Foldable side bars Safe bar
 System to lock the legs in closed position to lift and carry the device without using the wheels

Accessories

| | |
|----------|---|
| EN90003C | END-T Table for stretchers 10g cert. |
| ST00491A | STX 90 Telescopic headrest for stretchers |
| IF01044B | Track 4 telescopic IV holder for Ella Self |
| CB09025C | Tanker oxygen tank holder for stretcher side bar |
| CB09028C | Tanker support for Safe Bar |
| ST70000A | QMX777 Black PVC anatomic mattress |
| ST00497B | DNA STRAP Belt w/integrated re-winding system |
| ST00498B | DNA STRAP Thorax. Belt w/integrated re-winding system |
| ST42100A | FP Posterior fastener for self loading stretchers |
| ST42200A | FA Anterior fastener for self loading stretchers |

Standard Equipment

| | |
|----------|--------------------------------|
| ST00592A | STX 592 - TWO PIECE METAL BELT |
|----------|--------------------------------|

Technical data

| | |
|--------------------------------------|-------------------------------|
| Length | 1980 ±10 mm |
| Width | 570 ± 10 mm |
| Weight | 40 ± 1 kg |
| Loading capacity | 170 kg |
| Wheels | ø200 mm |
| Height of the loading carriage wheel | Min 44 ± 1cm Max 55 ± 1 cm |
| Frame material | Steel |
| Side bars length | 680 mm |
| Side bars Height | 200 mm |

Class I Med.Device according to Dir. 93/42/CEE

ELLA SELF T BLUE W/WHEELS Ø200

ST10600C



CND Classification V0899

Registration number 159744

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it