

ELLA SELF

Self loading stretcher



Spencer Ella Self is a compact and lightweight self-loading stretcher for transport of patients. It can be used for rescue operations on road and sanitary vehicles. The legs fold independently operating the specific commands and open automatically when unloading.

Specific Characteristics

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- Board with Trendelenburg antishock position
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bars Safe bar
- System to lock the legs in closed position to lift and carry the device without using the wheels
- Equipped with telescopic handles to make lifting operations easier

Accessories

EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
IF01044B	Track 4 telescopic IV holder for Ella Self
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar
ST70000A	QMX777 Black PVC anatomic mattress
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integrated re-winding system
ST42100A	FP Posterior fastener for self loading stretchers
ST42200A	FA Anterior fastener for self loading stretchers

Standard Equipment

ST00592A	STX 592 - TWO PIECE METAL BELT
----------	--------------------------------

Technical data

Length	1980 ±10 mm
Width	570 ± 10 mm
Weight	40 ± 1 kg
Loading capacity	170 kg
Wheels	ø200 mm
Height of the loading carriage wheel	Min 44 ± 1cm Max 55 ± 1 cm
Frame material	Steel
Side bars length	680 mm
Side bars Height	200 mm

Class I Med.Device according to Dir. 93/42/CEE

ELLA SELF T SILVER TELESCOPIC
HANDLES WHEELS °200M



ST10701C

CND Classification	V0899
Registration number	159730

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it