


NRT

Neil Robertson Type



Recovery stretcher designed for use in difficult conditions and in confined spaces.

Specific features

Low storage space required
 Transport handles allow easy use.
 Integrated adjustable belts for patient
 Belt for head stabilization

Technical data

Dimensions	1080x1830 ± 5 mm
Thickness	20 ± 5 mm
Material	Wood, Cotton
N° of handles	6
Loading capacity	220 kg
Weight	7 ± 0,5 kg

Class I Med.Device according to Dir. 93/42/CEE

NRT - NEIL ROBERTSON
 TYPERECOVERY STRETCHER
ST30150A



CND Classification	V9099
Identification number	216053

Rev.1 (15/09/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it