

Spencer 450

Transport chair with two wheels



Spencer 450 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific Characteristics

- Easy to clean PVC sheets
- 4 rear foldable handles to save space and to ensure lifting versatility
- Telescopic front handles to facilitate the handling
- Foldable to save space
- Strong and light frame

Reference standard

NA

Standard Equipment

ST30443A Pair of belts for chair 450/455

Accessories

ST30451A	MOUNTING BRACKETS FOR CHAIR SPENCER 450-55/420-25
ST36452E	UPPER FIXATION SYSTEM FOR TRANSPORT CHAIRS
ST40428B	Stabilizer stuffed seat/backrest profile series

Technical data

Length opened w/closed handles	610 ± 10 mm
Length opened w/opened handles	1140 ± 10 mm
Width	500 ± 10 mm
Height opened	920 ± 10 mm
Thickness folded	235 ± 10 mm
Rear wheels	ø125 mm
Materials	Steel/Al/PU/PVC
Weight	11,5 kg ± 0,3
Loading capacity	180 kg

Class I Med.Device according to Dir. 93/42/CEE

STAIR CHAIR SPENCER 450 WITH 2 WHEELS



ST30450A

CND Classification	V9099
Registration number	93460

Rev.3 (02/04/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
spencer.it