

Spencer 450/B

Transport chair with armrests



Spencer 450/B is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific Characteristics

Easy to clean PVC sheets
Foldable rear handles to reduce encumbrance.
Telescopic front handles to facilitate the handling
Slim when folded
Strong and light frame
Armrests to increase stability and comfort of the patient

Standard Equipment

ST30443A	Pair of belts for chair 450/ 455
----------	----------------------------------

Accessories

ST30451A	MOUNTING BRACKETS FOR CHAIR SPENCER 450-55/420-25
ST36452E	UPPER FIXATION SYSTEM FOR TRANSPORT CHAIRS
ST40428B	Stabilizer stuffed seat/backrest profile series

Technical data

Length opened w/closed handles	610 ± 10 mm
Length opened w/opened handles	1140 ± 10 mm
Width	510 ± 10 mm
Height opened	910 ± 10 mm
Height folded	920 ± 10 mm
Thickness folded	235 ± 10 mm
Rear wheels	ø125 mm
Materials	Acciaio/Al/PU/PVC
Weight	11,7 kg ± 0,3
Loading capacity	180 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 450/B TRANSPORT CHAIR
2WHEELS W/ARMREST



ST40450A

CND Classification V9099

Registration number 93461

Rev.2 (02/04/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it