

SPENCER 402-2

Compact transport chair with two wheels



Spencer 402 is the compact transport chair designed not only for those who have the problem of limited space in working space and in sanitary vehicle but also for those who want a versatile, dynamic device with both limited weight and dimensions.

The innovative kinematic opening and closing system is of fundamental importance for the extremely fast opening and closing of the device.

Extending front handles and rear handles with armrest make the chair ideal for patient transport on staircases.

The seat for the patient has been manufactured in robust PVC Spentex and is easily removed for washing and cleaning.

TECHNICAL DATA

| | |
|------------------------------|------------|
| Height opened | 860 mm |
| Height closed | 610 mm |
| Width | 520 mm |
| Depth folded | 160 mm |
| Length opened | 630 mm |
| Length with extended handles | 1020 mm |
| Surface treatment | Chromed |
| Materials | Steel, PVC |
| Weight | 9.2 kg |
| Loading capacity | 150 kg |
| Rear Wheels | Ø 125 mm |
| Latex free | |

ST00403B

SPENCER 402/2



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>