

## Spencer 406

Transport chair with 2 wheels



Spencer 406 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

### Specific features

Easy to clean PVC sheets
Robust and lightweight structure
Large rear wheels to reduce vibration during patient transfer
Quick and intuitive closure system
Safety system against accidental closure

### Reference standard

NA

### Standard equipment

ST00423A	Copule of belts for chairs
----------	----------------------------

### Optional accessories

ST42704E	Wall fastener 10g certified for chair 406
----------	---

### Technical data

Length	840 ± 10 mm
Width	410 ± 10 mm
Height opened	940 ± 10 mm
Height folded	770 ± 10 mm
Thickness folded	235 ± 10 mm
Rear wheels	∅ 200 mm
Materials	Steel/Al/PVC
Weight	10 kg ± 0,3
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

406 SPENCER AMBULANCE  
CARRYING CHAIR  
ST00406A



CND Classification	V9099
Registration number	93547

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)