

SPENCER 406B

Compact evacuation chair



The base model of the Compact Series offers high levels of safety and durability. The compactness, simplicity of use, rapidity of opening and closing movements, the solidity of the welded frame, the secure posture of the patient and the contained cost translate into dynamism, agility and capacity for immediate answers in any situation. Armrests and large posterior auxiliary handles increase comfort and maneuverability. The Spencer 406 chair is supplied with three polypropylene 50 mm [1.96 in] straps.

TECHNICAL DATA

Height opened	930 mm
Height folded	770 mm
Width	410 mm
Thickness folded	220 mm
Length	840 mm
Length with opened armrests	870 mm
Surface treatment	Chromed
Materials	Steel, PVC
Weight	11,6 kg
Load Capacity	150 kg
Rear Wheels	Ø 200x50 mm
Latex free	

SPENCER 406B



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it

Rev.1 (17/10/2013)

To verify the correctness of this document, please contact our Commercial Dept. at info@spencer.it